

Product Datasheet

Neisseria gonorrhoeae Antibody JP024

Product Name	Neisseria gonorrhoeae Antibody
Catalog Number	JP024
Species/Host	Rabbit
Tested Applications	LF, ELISA, WB, IHC, IFA
Immunogen	Neisseria gonorrhoeae whole cells
Brand	EastCoast Bio
Form/Appearance	Polyclonal Antibody, purified IgG Fraction
Storage	Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Note	Biotin, Fluorescein, and Peroxidase conjugates of this preparation are available.
Clonality	Polyclonal
Purity	>95%
Buffer	10 mM Phosphate Buffered Saline pH 7.2
Cross Reactivity	May react with related microbes (antiserum in unabsorbed)
Specificity	All Neisseria Gonorrhoeae Antigens
Preservatives	0.1% Sodium Azide