

Product Datasheet

Papillomavirus Antibody JP026

Product Name	Papillomavirus Antibody
Catalog Number	JP026
Species/Host	Goat
Tested Applications	LF, ELISA, WB, IHC, IFA
Immunogen	Bovine Papillomavirus
Brand	EastCoast Bio
Form/Appearance	Purified Polyclonal Antibody, IgG Fraction
Storage	Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Clonality	Polyclonal
Purity	>95%
Buffer	10 mM Phosphate Buffered Saline, pH 7.2
Cross Reactivity	Nonreactive with normal tissue and cultured cells
Specificity	Papillomavirus group including HPV types 11 & 16, Bovine Papillomavirus
Preservatives	0.1% Sodium Azide