

Product Datasheet

Parainfluenza Virus Antibody JP027

Product Name	Parainfluenza Virus Antibody
Catalog Number	JP027
Species/Host	Goat
Tested Applications	LF, ELISA, WB, IHC, IFA
Immunogen	Parainfluenza Virus, Type 3 (human isolate)
Brand	EastCoast Bio
Form/Appearance	Purified Polyclonal Antibody, IgG Fraction
Storage	Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Note	Biotin, Fluorescein, and Peroxidase conjugates of this preparation are available.
Clonality	Polyclonal
Purity	>95%
Buffer	10 mM Phosphate Buffered Saline, pH 7.2
Cross Reactivity	Cross-reacts with Parainfluenza Virus type 2, reactivity with Bovine type 3 and Canine Parainfluenza Virus, minimal reactivity with type 1, no reactivity with Hep-2 cells
Specificity	Structural proteins of Parainfluenza Virus, type 3
Preservatives	0.1% Sodium Azide