

Product Datasheet

Pichia pastoris Antibody JP028

Product Name	Pichia pastoris Antibody
Catalog Number	JP028
Species/Host	Rabbit
Tested Applications	LF, ELISA, WB, IHC, IFA
Immunogen	P. pastoris whole cells
Brand	EastCoast Bio
Form/Appearance	Purified Polyclonal Antibody, IgG Fraction
Storage	Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Clonality	Polyclonal
Purity	>95%
Buffer	10 mM Phosphate Buffered Saline, pH 7.2
Cross Reactivity	Strong reactivity with C. albicans and weak reactivity with S. cerevisiae
Specificity	All antigens of P. pastoris
Preservatives	0.1% Sodium Azide