

Product Datasheet

Poliovirus Antibody JP029

Product Name	Poliovirus Antibody
Catalog Number	JP029
Species/Host	Goat
Tested Applications	LF, ELISA, WB, IHC, IFA
Immunogen	Purified Poliovirus types 1-3
Brand	EastCoast Bio
Form/Appearance	Purified Polyclonal Antibody, IgG Fraction
Storage	Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Note	Biotin, Fluorescein, and Peroxidase conjugates of this preparation are available.
Clonality	Polyclonal
Purity	>95%
Buffer	10 mM Phosphate Buffered Saline, pH 7.2
Specificity	Poliovirus Types 1-3
Preservatives	0.1% Sodium Azide