

Product Datasheet

Saccharomyces cerevisiae Antibody JP031

Product Name	Saccharomyces cerevisiae Antibody
Catalog Number	JP031
Species/Host	Rabbit
Tested Applications	LF, ELISA, WB, IHC, IFA
Immunogen	S. cerevisiae whole cells
Brand	EastCoast Bio
Form/Appearance	Purified Polyclonal Antibody, IgG Fraction
Storage	Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Clonality	Polyclonal
Purity	>95%
Buffer	10 mM Phosphate Buffered Saline, pH 7.2
Cross Reactivity	Not determined
Specificity	All Antigens Of S. cerevisiae
Preservatives	0.1% Sodium Azide