

# Product Datasheet

## Streptococcus pneumoniae Antibody JP034

Product Name	Streptococcus pneumoniae Antibody
Catalog Number	JP034
Species/Host	Rabbit
Tested Applications	LF, ELISA, WB, IHC, IFA
Immunogen	Whole cell mix of multiple serotypes
Brand	EastCoast Bio
Form/Appearance	Purified Polyclonal Antibody, IgG Fraction
Storage	Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Note	Biotin, Fluorescein, and Peroxidase conjugates of this preparation are available.
Clonality	Polyclonal
Purity	>95%
Buffer	10 mM Phosphate Buffered Saline, pH 7.2
Cross Reactivity	May react with related microbes (antiserum in unabsorbed)
Specificity	Known Reactivity With S. Pneumoniae Types 3,4,6,7,9,14,18,19 & 23
Preservatives	0.1% Sodium Azide