

Product Datasheet

Tetanus toxoid Antibody JP035

Product Name	Tetanus toxoid Antibody
Catalog Number	JP035
Species/Host	Goat
Tested Applications	LF, ELISA, WB, IHC, IFA
Immunogen	Purified tetanus toxoid
Brand	EastCoast Bio
Form/Appearance	Purified Polyclonal Antibody, IgG Fraction
Storage	Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Clonality	Polyclonal
Purity	>95%
Buffer	10 mM Phosphate Buffered Saline pH 7.2
Specificity	Tetanus Toxoid
Preservatives	0.1% Sodium Azide