

# Product Datasheet

## Human Immunodeficiency Virus 1 (HIV 1-p24) Antibody-Biotin JP040

Product Name	Human Immunodeficiency Virus 1 (HIV 1-p24) Antibody-Biotin
Catalog Number	JP040
Species/Host	Rabbit
Tested Applications	ELISA
Immunogen	Native HIV antigen
Brand	EastCoast Bio
Form/Appearance	Purified Polyclonal Antibody-Biotin Conjugate
Storage	Short Term (≤ 2 weeks): 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Clonality	Polyclonal
Purity	>90%
Buffer	10 mM Phosphate Buffered Saline, pH 7.4 with 50% glycerol
Matched Pairs	HM083 or HM084