

Product Datasheet

Bovine Serum Albumin (BSA) Antibody JP041

Product Name	Bovine Serum Albumin (BSA) Antibody
Catalog Number	JP041
Description	Polyclonal Antibody, delipidated serum
Species/Host	Rabbit
Tested Applications	WB, ELISA, IP
Immunogen	Highly Purified Bovine Serum Albumin
Brand	EastCoast Bio
Form/Appearance	Polyclonal Antibody, delipidated serum
Storage	+2-8 °C. Do not freeze.
Clonality	Polyclonal
Specificity	Single precipitin arc using immunoelectrophoresis against BSA and the albumin in bovine serum
Preservatives	0.09% Sodium Azide