

Product Datasheet

Ferritin Antibody JP054

Product Name	Ferritin Antibody
Catalog Number	JP054
Description	Affinity purified and filtered complete IgY fraction
Species/Host	Chicken
Tested Applications	LF, ELISA
Immunogen	Ferritin isolated from human liver
Brand	EastCoast Bio
Form/Appearance	Affinity purified and filtered complete IgY fraction
Storage	+2-8 °C
Note	Filter prior to usage
Isotype	IgG
Clonality	Polyclonal
Buffer	0.1M Sodium phosphate buffer, pH 7.2 with 0.15M NaCl
Preservatives	0.09% Sodium Azide