

## Product Datasheet

### Ferritin Antibody JP054

<b>Product Name</b>	Ferritin Antibody
<b>Catalog Number</b>	JP054
<b>Description</b>	Affinity purified and filtered complete IgY fraction
<b>Species/Host</b>	Chicken
<b>Tested Applications</b>	ELISA, LF
<b>Immunogen</b>	Ferritin isolated from human liver
<b>Brand</b>	EastCoast Bio
<b>Form/Appearance</b>	Affinity purified and filtered complete IgY fraction
<b>Storage</b>	+2-8 °C
<b>Note</b>	Filter prior to usage
<b>Isotype</b>	IgG
<b>Clonality</b>	Polyclonal
<b>Buffer</b>	0.1M Sodium phosphate buffer, pH 7.2 with 0.15M NaCl
<b>Preservatives</b>	0.09% Sodium Azide