

Product Datasheet

C-Reactive Protein, Canine (cCRP) Antibody JP057

Product Name	C-Reactive Protein, Canine (cCRP) Antibody
Catalog Number	JP057
Description	Purified complete IgY fraction
Species/Host	Chicken
Tested Applications	ELISA, LF
Immunogen	Canine CRP isolated from dog serum
Brand	EastCoast Bio
Form/Appearance	Purified complete IgY fraction
Storage	+2-8 °C
Clonality	Polyclonal
Buffer	0.1M Sodium Phosphate Buffer, pH 7.2 with 150 mM Sodium Chloride
Preservatives	0.09% Sodium Azide