

Product Datasheet

Human CD54 ELISA Set 851.590

| | |
|---------------------|--|
| Product Name | Human CD54 ELISA Set |
| Catalog Number | 851.590 |
| Description | Diaclone ELISA Sets include capture and biotinylated detection antibody, one standard per plate, Streptavidin-HRP, TMB, and detailed procedure including buffer composition. |
| Range | 0.25-8 ng/ml |
| Sensitivity | 0.1 ng/ml |
| Tested Applications | ELISA |
| Alternative Names | ICAM-1 |
| Brand | Diaclone |
| Storage | Store at +4°C |
| Cross Reactivity | No cross reactivity with other human adhesion molecules |
| Specificity | Recognizes both natural and recombinant human soluble ICAM-1 |
| Sample Types | Serum, Plasma, Cell culture supernatant |