

# Product Datasheet

## Human CD54 ELISA Set 851.590.001

Product Name	Human CD54 ELISA Set
Catalog Number	851.590.001
Description	Diaclone ELISA Sets include capture and biotinylated detection antibody, one standard per plate, Streptavidin-HRP, TMB, and detailed procedure including buffer composition.
Range	0.25-8 ng/ml
Sensitivity	0.1 ng/ml
Tested Applications	ELISA
Alternative Names	ICAM-1
Brand	Diaclone
Storage	Store at +4°C
Cross Reactivity	No cross reactivity with other human adhesion molecules
Specificity	Recognizes both natural and recombinant human soluble ICAM-1
Sample Types	Serum, Plasma, Cell culture supernatant