

## Product Datasheet

### Apolipoprotein A1 (ApoA1) Antibody JP067

<b>Product Name</b>	Apolipoprotein A1 (ApoA1) Antibody
<b>Catalog Number</b>	JP067
<b>Species/Host</b>	Goat
<b>Tested Applications</b>	TIA
<b>Immunogen</b>	Purified Human Apolipoprotein A1
<b>Brand</b>	EastCoast Bio
<b>Form/Appearance</b>	Polyclonal, Delipidated and 0.2 µm filtered
<b>Storage</b>	Short Term (≤ 2 weeks): 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
<b>Note</b>	Treated for >2 hours @ pH < 5.0
<b>Clonality</b>	Polyclonal
<b>Buffer</b>	20 mM Tris at pH 7.4 with 250mM NaCl
<b>Specificity</b>	Single precipitin arc visible with immunoelectrophoresis
<b>Preservatives</b>	0.1% Sodium Azide