

Product Datasheet

Apolipoprotein B (ApoB) Antibody JP068

Product Name	Apolipoprotein B (ApoB) Antibody
Catalog Number	JP068
Species/Host	Goat
Tested Applications	TIA
Immunogen	Purified Human Apolipoprotein B
Brand	EastCoast Bio
Form/Appearance	Polyclonal, Delipidated and 0.2 µm filtered
Storage	Short Term (≤ 2 weeks): 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Note	Treated for >2 hours @ pH < 5.0
Clonality	Polyclonal
Buffer	20 mM Tris, pH 7.4 with 250 mM NaCl
Specificity	Single precipitin arc visible with immunoelectrophoresis
Preservatives	0.1% Sodium Azide