

Product Datasheet

Transferrin Antibody JP072

Product Name	Transferrin Antibody
Catalog Number	JP072
Species/Host	Goat
Tested Applications	TIA
Immunogen	Purified Human Transferrin
Brand	EastCoast Bio
Form/Appearance	Polyclonal, Delipidated and 0.2 µm filtered
Storage	Short Term (≤ 2 weeks): 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Note	Treated for > 2 hours @ pH <5.0
Clonality	Polyclonal
Buffer	20 mM Tris, pH 7.4 with 250mM NaCl
Specificity	Single precipitin arc visible with immunoelectrophoresis
Preservatives	0.1% Sodium Azide