

Product Datasheet

Transferrin Antibody JP072

| | |
|----------------------------|---|
| Product Name | Transferrin Antibody |
| Catalog Number | JP072 |
| Species/Host | Goat |
| Tested Applications | TIA |
| Immunogen | Purified Human Transferrin |
| Brand | EastCoast Bio |
| Form/Appearance | Polyclonal, Delipidated and 0.2 µm filtered |
| Storage | Short Term (≤ 2 weeks): 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing. |
| Note | Treated for > 2 hours @ pH <5.0 |
| Clonality | Polyclonal |
| Buffer | 20 mM Tris, pH 7.4 with 250mM NaCl |
| Specificity | Single precipitin arc visible with immunoelectrophoresis |
| Preservatives | 0.1% Sodium Azide |