

Product Datasheet

Apolipoprotein Lp(a) Antibody JP078

Product Name	Apolipoprotein Lp(a) Antibody
Catalog Number	JP078
Description	Polyclonal, Delipidated and 0.2 µm filtered
Species/Host	Goat
Tested Applications	TIA
Immunogen	Human Apolipoprotein (a)
Brand	EastCoast Bio
Form/Appearance	Polyclonal, Delipidated and 0.2 µm filtered
Storage	Short Term (≤ 2 weeks): +2-8 °C. Long Term: -20 °C. Avoid repeated freezing and thawing.
Note	Filtration prior to use is recommended if particulates are observed
Clonality	Polyclonal
Buffer	0.05M Tris, pH 7.5 with 0.5M NaCl
Specificity	Single precipitin arc visible with immunoelectrophoresis
Preservatives	0.09 % Sodium Azide