

Product Datasheet

Coronavirus, Canine Antibody JP082

Product Name	Coronavirus, Canine Antibody
Catalog Number	JP082
Description	Affinity Purified Polyclonal Antibody
Species/Host	Goat
Tested Applications	WB, ELISA
Immunogen	Canine Coronavirus nucleocapsid protein
Brand	EastCoast Bio
Form/Appearance	Affinity Purified Polyclonal Antibody
Storage	+2-8 °C. Do not freeze.
Clonality	Polyclonal
Buffer	10 mM Phosphate Buffered Saline, pH 7.2
Preservatives	0.1% Sodium Azide