

## Product Datasheet

### Streptococcus Group A Antibody JP086

|                            |  |
|----------------------------|--|
| <b>Product Name</b>        | Streptococcus Group A Antibody   |
| <b>Catalog Number</b>      | JP086  |
| <b>Description</b>         | Affinity Purified Polyclonal Antibody  |
| <b>Species/Host</b>        | Rabbit   |
| <b>Tested Applications</b> | ELISA, LF  |
| <b>Immunogen</b>           | S. pyogenes  |
| <b>Brand</b>               | EastCoast Bio  |
| <b>Form/Appearance</b>     | Affinity Purified Polyclonal Antibody  |
| <b>Storage</b>             | +2-8 °C  |
| <b>Clonality</b>           | Polyclonal   |
| <b>Purity</b>              | ≥ 95 %   |
| <b>Buffer</b>              | 0.02M Sodium Phosphate, 0.15M Sodium Chloride  |
| <b>Cross Reactivity</b>    | No reactivity with tested strains of Streptococcus Group B, C, F, G, S. aureus, H. influenzae, N. perflava |
| <b>Specificity</b>         | Group A Streptococcus  |
| <b>Preservatives</b>       | 0.05% Sodium Azide   |