

## Product Datasheet

### Cystatin C Antibody JP089

<b>Product Name</b>	Cystatin C Antibody
<b>Catalog Number</b>	JP089
<b>Species/Host</b>	Rabbit
<b>Tested Applications</b>	ELISA
<b>Brand</b>	EastCoast Bio
<b>Form/Appearance</b>	Liquid; Purified Polyclonal Antibody
<b>Storage</b>	Short Term ( $\leq 2$ weeks): 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
<b>Clonality</b>	Polyclonal
<b>Purity</b>	>90%
<b>Buffer</b>	10 mM Phosphate Buffered Saline, pH 7.4