
Product Datasheet

Cystatin C Antibody JP089

Product Name	Cystatin C Antibody
Catalog Number	JP089
Species/Host	Rabbit
Tested Applications	ELISA
Brand	EastCoast Bio
Form/Appearance	Liquid; Purified Polyclonal Antibody
Storage	Short Term (≤ 2 weeks): 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Clonality	Polyclonal
Purity	>90%
Buffer	10 mM Phosphate Buffered Saline, pH 7.4