

Product Datasheet

Beta-2 Microglobulin (B2M) Antibody JP090

Product Name	Beta-2 Microglobulin (B2M) Antibody
Catalog Number	JP090
Species/Host	Rabbit
Tested Applications	ELISA
Immunogen	Human Beta-2 Microglobulin
Brand	EastCoast Bio
Form/Appearance	Affinity Purified Polyclonal Antibody
Storage	Short Term (≤ 2 weeks): 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Clonality	Polyclonal
Purity	>90%
Buffer	10 mM Phosphate Buffered Saline, pH 7.2
Preservatives	0.09% Sodium Azide