

Product Datasheet

6-Monoacetylmorphine (6-MAM) Antibody JP091

Product Name	6-Monoacetylmorphine (6-MAM) Antibody
Catalog Number	JP091
Description	Protein A Purified Polyclonal Antibody, Ig Fraction
Species/Host	Rabbit
Tested Applications	ELISA, LF
Immunogen	6-Monoacetylmorphine-BgG
Alternative Names	Heroin
Brand	EastCoast Bio
Form/Appearance	Protein A Purified Polyclonal Antibody, Ig Fraction
Storage	Short Term (≤ 2 weeks): +2-8 °C. Long Term: -20 °C. Avoid repeated freezing and thawing.
Clonality	Polyclonal
Buffer	0.05M Phosphate, pH 7.15 with 0.85% Sodium Chloride
Cross Reactivity	Desomorphine <0.1%, Diacetylmorphine <2.5%
Specificity	6-Monoacetylmorphine
Matched Pairs	170031 (capture)
Preservatives	0.05% Sodium Azide