

## Product Datasheet

### 6-Monoacetylmorphine (6-MAM) Antibody JP092

<b>Product Name</b>	6-Monoacetylmorphine (6-MAM) Antibody
<b>Catalog Number</b>	JP092
<b>Description</b>	Protein A Purified Polyclonal Antibody, Ig Fraction
<b>Species/Host</b>	Rabbit
<b>Tested Applications</b>	ELISA, LF
<b>Immunogen</b>	6-Monoacetylmorphine-BgG
<b>Alternative Names</b>	Heroin
<b>Brand</b>	EastCoast Bio
<b>Form/Appearance</b>	Protein A Purified Polyclonal Antibody, Ig Fraction
<b>Storage</b>	Short Term ( $\leq 2$ weeks): +2-8 °C. Long Term: -20 °C. Avoid repeated freezing and thawing.
<b>Clonality</b>	Polyclonal
<b>Buffer</b>	0.05M Phosphate, pH 7.15 with 0.85% Sodium Chloride
<b>Cross Reactivity</b>	Desomorphine <0.1%, Diacetylmorphine <2.5%, Morphine <0.1%, Nalorphine <0.1%
<b>Specificity</b>	6-Monoacetylmorphine
<b>Matched Pairs</b>	170031 (capture)
<b>Preservatives</b>	0.05% Sodium Azide