

Product Datasheet

Human Immunoglobulin E (IgE) Antibody JP094

Product Name	Human Immunoglobulin E (IgE) Antibody
Catalog Number	JP094
Description	Affinity Purified Polyclonal Antibody
Species/Host	Goat
Tested Applications	WB, IHC, ELISA
Immunogen	Purified Human IgE Heavy (e) chain
Brand	EastCoast Bio
Form/Appearance	Affinity Purified Polyclonal Antibody
Storage	+2-8 °C
Clonality	Polyclonal
Purity	≥ 95 % by SDS-PAGE
Buffer	10 mM Sodium Phosphate, pH 7.2 with 0.15 M Sodium Chloride
Cross Reactivity	No reactivity with Human immunoglobulin light chains or other Human serum proteins
Specificity	Human IgE Heavy (e) chain
Preservatives	0.05% Sodium Azide