

Product Datasheet

Human Kappa Light Chain Antibody-AP JP095

Product Name	Human Kappa Light Chain Antibody-AP
Catalog Number	JP095
Species/Host	Goat
Tested Applications	ELISA, WB, IHC, IF
Immunogen	Purified Human Kappa (κ) Light chain
Brand	EastCoast Bio
Form/Appearance	Affinity Purified Polyclonal Antibody
Storage	2-8°C
Clonality	Polyclonal
Purity	$\geq 95\%$ by SDS-PAGE
Buffer	10 mM Sodium Phosphate, pH 7.2 with 0.15 M Sodium Chloride
Cross Reactivity	No reactivity with Human immunoglobulin heavy chains, Human immunoglobulin lambda (λ) light chains or other Human serum proteins
Specificity	Human Kappa (κ) Light chains
Preservatives	0.05% Sodium Azide