

Product Datasheet

K2 Antibody JP098

Product Name	K2 Antibody
Catalog Number	JP098
Description	Protein A Purified Polyclonal Antibody
Species/Host	Rabbit
Tested Applications	ELISA, LF
Immunogen	K2-BgG
Brand	EastCoast Bio
Form/Appearance	Protein A Purified Polyclonal Antibody
Storage	Short Term: +2-8 °C. Long Term: -20 °C. Avoid repeated freezing and thawing.
Clonality	Polyclonal
Buffer	0.05 M Phosphate Buffer, pH 7.15 with 0.85% Sodium Chloride
Cross Reactivity	JWH-018 100%; AM-2232 35.7%; JWH-073 10.1%; JWH-200 3.7%; JWH-019 3.3; JWH-203 <0.2%, JWH-147, JWH-175, HU-210, CP-47,947, AH-7921 <0.1%
Specificity	K2
Matched Pairs	LA288 (capture)
Preservatives	0.05% Sodium Azide