

Product Datasheet

Hydrocodone Antibody JP101

Product Name	Hydrocodone Antibody
Catalog Number	JP101
Description	Protein A Purified Polyclonal Antibody
Species/Host	Rabbit
Tested Applications	ELISA, LF
Immunogen	Hydrocodone-BgG
Brand	EastCoast Bio
Form/Appearance	Protein A Purified Polyclonal Antibody
Storage	Short Term: +2-8 °C. Long Term: -20 °C. Avoid repeated freezing and thawing.
Clonality	Polyclonal
Buffer	0.05 M Phosphate Buffer, pH 7.15 with 0.85% Sodium Chloride
Cross Reactivity	Hydrocodone 100%; Hydromorphone 2.2%; Codeine <0.2%; Morphine, Oxycodone <0.1%
Specificity	Hydrocodone
Matched Pairs	LA289 (capture)
Preservatives	0.05% Sodium Azide