

# Product Datasheet

## C-Reactive Protein (CRP) Antibody JP107

|                     |                                                          |
|---------------------|----------------------------------------------------------|
| Product Name        | C-Reactive Protein (CRP) Antibody                        |
| Catalog Number      | JP107                                                    |
| Species/Host        | Goat                                                     |
| Tested Applications | ELISA, IEP, TIA                                          |
| Brand               | EastCoast Bio                                            |
| Form/Appearance     | Polyclonal, Delipidated and 0.2 µm filtered              |
| Storage             | 2-8°C                                                    |
| Clonality           | Polyclonal                                               |
| Buffer              | 0.2M Tris, 0.5M NaCl, pH 7.2                             |
| Specificity         | Single precipitin arc visible with immunoelectrophoresis |
| Preservatives       | 0.05% Sodium Azide                                       |