

Product Datasheet

Human Immunoglobulin M (IgM) Scavenger Antibody JP111

Product Name	Human Immunoglobulin M (IgM) Scavenger Antibody
Catalog Number	JP111
Species/Host	Mouse
Immunogen	Human IgM (non-specific)
Brand	EastCoast Bio
Form/Appearance	Purified Monoclonal Antibody
Storage	Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Clonality	Monoclonal
Buffer	0.015M Potassium Phosphate, pH 7.2 with 0.85% Sodium Chloride
Preservatives	0.05% Sodium Azide