

# Product Datasheet

## Bacillus Antibody, Biotin JP112

Product Name	Bacillus Antibody, Biotin
Catalog Number	JP112
Species/Host	Rabbit
Tested Applications	ELISA, WB, IHC, IFA
Immunogen	Purified spores of B. cereus and B. subtilis
Brand	EastCoast Bio
Form/Appearance	Purified Polyclonal Antibody, Biotin Conjugated
Storage	Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Clonality	Polyclonal
Purity	>95%
Buffer	10 mM Phosphate Buffered Saline, pH 7.2
Cross Reactivity	May react with multiple Bacillus species (antiserum is unabsorbed)
Preservatives	0.1% Sodium Azide