

Product Datasheet

Mycobacterium tuberculosis (TB) Antibody, FITC JP123

Product Name	Mycobacterium tuberculosis (TB) Antibody, FITC
Catalog Number	JP123
Species/Host	Rabbit
Tested Applications	IFA
Immunogen	Purified protein derivative (PPD) of M. tuberculosis
Brand	EastCoast Bio
Form/Appearance	Purified Polyclonal Antibody, Fluorescein Conjugate
Storage	Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Note	Protect from light
Isotype	IgG
Clonality	Polyclonal
Purity	>95%
Buffer	10 mM Phosphate Buffered Saline pH 7.2 with 10 mg/ml BSA
Cross Reactivity	Reactive with several related microbes including M. avium, M. phlei, and M. parafortuitum (antiserum is unabsorbed)
Specificity	M. tuberculosis
Preservatives	0.1% Sodium Azide