

Product Datasheet

Complement Component 3 (C3) Antibody JP127

Product Name	Complement Component 3 (C3) Antibody
Catalog Number	JP127
Species/Host	Goat
Tested Applications	TIA
Immunogen	Native Human C3
Brand	EastCoast Bio
Form/Appearance	Polyclonal, Delipidated and 0.2 µm filtered
Storage	2-8°C
Clonality	Polyclonal
Buffer	20 mM Tris, pH 7.4 with 500 mM NaCl
Specificity	Single precipitin arc visible with immunoelectrophoresis
Preservatives	0.1% Sodium Azide