Medix Biochemica

Oy Medix Biochemica Ab

Headquarters Klovinpellontie 3 FI-02180 Espoo

Finland

Email: sales@medixbiochemica.com

Medix Biochemica USA Inc.

10850 Metro Court Maryland Heights, MO 63043 United States of America

Email: sales@medixbiochemica.com

Product Datasheet

Histone H1 Antibody JP130

Product

Histone H1 Antibody

Name

Catalog Number JP130

Description

Purified Polyclonal Antibody

Species/Host Sheep

Tested

WB, IHC, IP

Applications

Immunogen

Bovine histone H1-RNA complex

Brand

EastCoast Bio

Form/AppearaPurified Polyclonal Antibody

Storage

Short Term (≤ 1 weeks): +2-8 °C. Long Term: -20 °C. Avoid

repeated freezing and thawing.

Clonality

Polyclonal

Buffer

10 mM Phosphate Buffered Saline

Cross

Reactive with many species

Reactivity