

# Product Datasheet

## Histone H1 Antibody JP130

Product Name	Histone H1 Antibody
Catalog Number	JP130
Species/Host	Sheep
Tested Applications	IP, WB, IHC
Immunogen	Bovine histone H1-RNA complex
Brand	EastCoast Bio
Form/Appearance	Purified Polyclonal Antibody
Storage	Short Term ( $\leq$ 1 weeks): 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Clonality	Polyclonal
Buffer	10 mM Phosphate Buffered Saline
Cross Reactivity	Reactive with many species