

Product Datasheet

Histone H2B Antibody JP131

Product Name	Histone H2B Antibody
Catalog Number	JP131
Species/Host	Sheep
Tested Applications	IP, WB
Immunogen	Bovine histone H2B-RNA complex
Brand	EastCoast Bio
Form/Appearance	Purified Polyclonal Antibody
Storage	Short Term (≤ 1 weeks): 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Clonality	Polyclonal
Buffer	10 mM Phosphate Buffered Saline
Cross Reactivity	Reactive with many species