

## Product Datasheet

### Histone H3 Antibody JP132

<b>Product Name</b>	Histone H3 Antibody
<b>Catalog Number</b>	JP132
<b>Species/Host</b>	Sheep
<b>Tested Applications</b>	IP, WB, IHC
<b>Immunogen</b>	Bovine histone H3-RNA complex
<b>Brand</b>	EastCoast Bio
<b>Form/Appearance</b>	Purified Polyclonal Antibody
<b>Storage</b>	Short Term ( $\leq 1$ weeks): 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
<b>Clonality</b>	Polyclonal
<b>Buffer</b>	10 mM Phosphate Buffered Saline
<b>Cross Reactivity</b>	Reactive with many species