Medix Biochemica

Oy Medix Biochemica Ab

Headquarters Klovinpellontie 3 Fl-02180 Espoo

Finland

Email: sales@medixbiochemica.com

Medix Biochemica USA Inc.

10850 Metro Court Maryland Heights, MO 63043 United States of America

Email: sales@medixbiochemica.com

Product Datasheet

Histone H3 Antibody JP132

Product Name Histone H3 Antibody

JP132

Catalog Number

Description Purified Polyclonal Antibody

Species/Host Sheep

Tested

WB, IHC, IP

Applications

Immunogen

Bovine histone H3-RNA complex

Brand

EastCoast Bio

Form/AppearaPurified Polyclonal Antibody

Short Term (≤ 1 weeks): +2-8 °C. Long Term: -20 °C. Avoid

repeated freezing and thawing.

Clonality

Storage

Polyclonal

Buffer

10 mM Phosphate Buffered Saline

Cross

Reactive with many species

Reactivity