

Product Datasheet

Streptococcus Group B Antibody JP133

Product Name	Streptococcus Group B Antibody
Catalog Number	JP133
Species/Host	Rabbit
Tested Applications	LF, ELISA, WB, IHC, IFA
Immunogen	Group B streptococci
Brand	EastCoast Bio
Form/Appearance	Purified Polyclonal Antibody, IgG Fraction
Storage	Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Note	Biotin, Fluorescein, and Peroxidase conjugates of this antiserum preparation are available.
Clonality	Polyclonal
Purity	>95%
Buffer	10 mM Phosphate Buffered Saline, pH 7.2
Cross Reactivity	Does not cross-react with other Streptococcus groups
Specificity	Group B Streptococci
Matched Pairs	JP133 (self-pairs)
Preservatives	0.1% Sodium Azide