

## Product Datasheet

### Enterovirus Antibody JP134

<b>Product Name</b>	Enterovirus Antibody
<b>Catalog Number</b>	JP134
<b>Description</b>	Purified Polyclonal Antibody
<b>Species/Host</b>	Rabbit
<b>Tested Applications</b>	ELISA, LF
<b>Brand</b>	EastCoast Bio
<b>Form/Appearance</b>	Purified Polyclonal Antibody
<b>Storage</b>	+2-8 °C
<b>Clonality</b>	Polyclonal
<b>Purity</b>	> 95 % by SDS-PAGE
<b>Buffer</b>	10 mM Phosphate Buffered Saline, pH 7.4
<b>Matched Pairs</b>	JP134 (self-pairs)
<b>Preservatives</b>	0.09% Sodium Azide