

Product Datasheet

S100A12 Antibody JP135

Product Name	S100A12 Antibody
Catalog Number	JP135
Description	Purified Polyclonal Antibody
Species/Host	Rabbit
Tested Applications	ELISA, LF
Brand	EastCoast Bio
Form/Appearance	Purified Polyclonal Antibody
Storage	+2-8 °C
Clonality	Polyclonal
Purity	> 95 % by SDS-PAGE
Buffer	10 mM Phosphate Buffered Saline, pH 7.4
Matched Pairs	JP136 (capture or detection)
Preservatives	0.09% Sodium Azide