

## Product Datasheet

### Elastase Antibody JP137

<b>Product Name</b>	Elastase Antibody
<b>Catalog Number</b>	JP137
<b>Description</b>	Affinity Purified Polyclonal Antibody
<b>Species/Host</b>	Rabbit
<b>Tested Applications</b>	ELISA, LF
<b>Brand</b>	EastCoast Bio
<b>Form/Appearance</b>	Affinity Purified Polyclonal Antibody
<b>Storage</b>	Short Term: +2-8 °C. Long Term: -20 °C. Avoid repeated freezing and thawing.
<b>Clonality</b>	Polyclonal
<b>Purity</b>	≥ 95 %
<b>Buffer</b>	Phosphate Buffered Saline, pH 7.4
<b>Matched Pairs</b>	JP138 (capture or detection)
<b>Preservatives</b>	0.09% Sodium Azide