

Product Datasheet

Elastase Antibody JP138

Product Name	Elastase Antibody
Catalog Number	JP138
Description	Protein A Purified Polyclonal Antibody
Species/Host	Rabbit
Tested Applications	ELISA, LF
Brand	EastCoast Bio
Form/Appearance	Protein A Purified Polyclonal Antibody
Storage	+2-8 °C
Clonality	Polyclonal
Purity	> 95 % by SDS-PAGE
Buffer	10 mM Phosphate Buffered Saline, pH 7.4
Matched Pairs	JP137 (capture or detection)
Preservatives	0.09% Sodium Azide