

Product Datasheet

Human Kappa Light Chain Antibody JP140

Product Name	Human Kappa Light Chain Antibody
Catalog Number	JP140
Description	Polyclonal, Delipidated and 0.2 µm filtered
Species/Host	Goat
Tested Applications	ELISA
Immunogen	Purified Human Kappa Light chain
Brand	EastCoast Bio
Form/Appearance	Polyclonal, Delipidated and 0.2 µm filtered
Storage	Short Term: +2-8 °C. Long Term: -20 °C. Avoid repeated freezing and thawing.
Clonality	Polyclonal
Buffer	TBS, pH 7.4
Specificity	Single precipitin arc visible with immunoelectrophoresis
Preservatives	0.1% Sodium Azide