

Product Datasheet

Pseudomonas Antibody, FITC JP143

Product Name	Pseudomonas Antibody, FITC
Catalog Number	JP143
Species/Host	Rabbit
Tested Applications	ELISA, WB, IHC, IFA
Immunogen	P. aeruginosa whole cells
Brand	EastCoast Bio
Form/Appearance	Polyclonal Antibody-FITC conjugate
Storage	Short Term: 2-8°C. Long Term: -20°C. Avoid repeated freezing and thawing.
Clonality	Polyclonal
Purity	>95%
Buffer	10 mM Phosphate Buffered Saline pH 7.2
Cross Reactivity	May react with related microbes (antiserum in unabsorbed)
Specificity	P. aeruginosa Soluble And Structural Proteins
Preservatives	0.1% Sodium Azide