

Product Datasheet

Prealbumin Antibody JP144

Product Name	Prealbumin Antibody
Catalog Number	JP144
Description	Polyclonal Antibody, defibrinated, delipidated, absorbed by solid phase chromatography as required and 0.2 µm filtered
Species/Host	Goat
Tested Applications	ELISA
Immunogen	Purified Human Prealbumin
Brand	EastCoast Bio
Form/Appearance	Polyclonal Antibody, defibrinated, delipidated, absorbed by solid phase chromatography as required and 0.2 µm filtered
Storage	Short Term (≤ 2 weeks): +2-8 °C. Long Term: -20 °C. Avoid repeated freezing and thawing.
Clonality	Polyclonal
Buffer	0.05M Tris, pH 7.5 with 0.5M NaCl
Specificity	Human
Preservatives	0.09% Sodium Azide