

## Product Datasheet

### Estradiol (E2) Antibody JP145

<b>Product Name</b>	Estradiol (E2) Antibody
<b>Catalog Number</b>	JP145
<b>Description</b>	Purified Polyclonal Antibody, IgG Fraction
<b>Species/Host</b>	Rabbit
<b>Tested Applications</b>	WB, IHC, ELISA
<b>Brand</b>	EastCoast Bio
<b>Form/Appearance</b>	Purified Polyclonal Antibody, IgG Fraction
<b>Storage</b>	+2-8 °C
<b>Clonality</b>	Polyclonal
<b>Buffer</b>	50 mM Sodium Phosphate Buffer with 0.9% NaCl, pH 7.3
<b>Cross Reactivity</b>	17 $\beta$ Estradiol 100%; Estrone <1%; Testosterone, Corticosterone, Pregnenolone, Androstenedione, Progesterone, Estriol < 0.1%
<b>Specificity</b>	Estradiol
<b>Preservatives</b>	0.025% Proclin 300